

# Dialysis Medical Solutions, LLC.

## Acid Concentrates

### 35X Dilution Acid Concentrate.

Product Code	Na+ (mEq/L)	Ca++ (mEq/L)	K+ (mEq/L)	Mg++ (mEq/L)	Cl- (mEq/L)	HOAc (mEq/L)	Dex (mg/dL)
AC-105	103.0	2.50	1.00	1.00	107.50	3.00	200.0
AC-106	105.0	2.50	2.00	1.00	110.50	3.00	200.0

### 36.83X Dilution Acid Concentrate.

Product Code	Na+ (mEq/L)	Ca++ (mEq/L)	K+ (mEq/L)	Mg++ (mEq/L)	Cl- (mEq/L)	HOAc (mEq/L)	Dex (mg/dL)
AC-100	81.0	2.00	2.50	1.00	86.50	4.00	100.0
AC-101	80.0	2.50	2.00	1.00	85.5	4.00	200.0
AC-102	80.0	3.00	2.00	1.00	86.0	4.00	200.0

### 45X Dilution Acid Concentrate.

Product Code	Na+ (mEq/L)	Ca++ (mEq/L)	K+ (mEq/L)	Mg++ (mEq/L)	Cl- (mEq/L)	HOAc (mEq/L)	Dex (mg/dL)
AC-103	100.0	2.50	2.00	0.75	105.25	4.00	200.0
AC-104	100.0	3.00	2.00	0.75	105.75	4.00	200.0